




Drill Pin Drivers

Designed to seamlessly attach to a powered driver, the Enztec Drill Pin Drivers are the ideal tool for efficient and secure pin installation.

The unique live spring mechanism, ensures reliable retention of tri-flat pins across a range of applications, providing confidence with every use.



Features

- Light weight and easy to use
- Ensures secure pin grip
- Easy attachment to powered tools (Zimmer or AO)
- Pull-to-release feature for convenient removal from the pin
- Compatible with tri flat pins

Specifications

- Materials: Stainless Steel
- Zimmer – 75mm long, 20g
- AO – 98mm long, 33g