


Partial Knee Impactor

Providing femoral trial extraction as well as impaction of final components.

The Enztec Partial Knee Impactor builds on the proven design of the Enztec Total Knee Femoral Impactor—one of our most recognised instruments, with over 15 years of clinical history.

Designed specifically for PKA/UKA, it features a polymer impaction head with three points of contact giving surgeons greater control and stability during guided placement and final impaction.

Its open design and removable impaction head ensure easy cleaning, delivering reliable performance, durability, and longevity.



Features

- Intuitive to use
- Configurable to suit your femoral components
- Screwable impaction handle allows wide range of femoral components to be firmly gripped
- Open design ensures it is easy to clean
- Replaceable impaction head

Specifications

- Materials: Stainless steel, Propylux
- Weight: 365g / 0.8lbs
- Length: 215mm / 8.4in
- Add your logo: available