

Straight Cup Impactor


Maximize in-situ visualisation.

Lightweight and well-balanced, the Enztec Straight Cup Impactor features a quick-connect trinket for easy detachment from the cup, allowing visual inspection or intraoperative imaging to verify cup position before final impaction.

Its simple, intuitive design is compatible with a wide range of acetabular cups or can be adapted to fit specific cups, ensuring seamless integration. Trinkets are interchangeable with the Enztec Classic Offset Cup Impactor X. The impactor also includes an alignment array connector that is compatible with Enztec's Alignment Arrays and NAVBIT SPRINT™.

Trust the Enztec Straight Cup Impactor to maximize cup placement control and enhance visualisation.





COMPATIBLE IMPLANTS

- Smith + Nephew R3 and REDAPT
- DePuy Pinnacle and Emphasys
- Stryker Trident
- ZB G7
- Signature Logical G
- Microport Procotyl P
- Corin Trinity
- + more

Features

- Allows in-situ visualization and intraoperative imaging of the cup before final impaction
- Compatible with Enztec's Alignment Arrays and NAVBIT SPRINT™
- Trinkets interchangeable with the Enztec Classic Offset Cup Impactor X

Specifications

- Materials: Stainless Steel, Propylux, Silicone
- Weight: 720g / 1.6lbs
- Length: 410mm / 16.1 in
- Add your logo: available.